

UltraComfort® - Asepsis 21® Operatory



We have tried our very best to accurately represent product colors on this brochure. However, due to the nature of printed media, the color of actual products may differ slightly from how they appear on your particular printed piece. In order to confirm your color selection, we strongly recommend that you contact Midmark Customer Service at 1-800-MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 Certified Company. For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495 Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214 or visit our website at midmark.com

© 2009 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0802-00 Rev. A1 (8/09)



A S E P S I S 2 1®
D E L I V E R Y S Y S T E M



R E L I A B L E



Midmark's patented Kink Valve technology offers an independent control system for each handpiece and an unprecedented level of reliability and control. Unlike traditional valve blocks that require ongoing maintenance of O rings and diaphragms, this system requires no maintenance.

The Asepsis 21[®] delivery system provides quality and reliability in a sleek, compact, value-priced package. Designed for improved asepsis, the air and water controls are protected from contamination beneath an easily removable top. Surfaces are smooth and easy to clean. Available in cart, console mount or left/right delivery, the Asepsis delivery unit will complement the way you work and deliver exceptional reliability.

Asepsis 21[®] Delivery Carts




Asepsis 21[®] Left/Right Duo Cart

- Three automatic Kink Valve handpiece positions
- Assistant's instrumentation: HVE, saliva ejector and syringe
- Water quick connect
- 21" x 24" stain resistant work surface
- Choice of 9" x 13" tray or tub support
- Adjustable height "H" frame (epoxy powder coat finish)



Asepsis 21[®] Assistant's Cart

- Assistant's instrumentation: HVE, saliva ejector and syringe
- Water quick connect
- 21" x 24" stain resistant work surface
- Choice of 9" x 13" tray or tub support
- Adjustable height "H" frame (epoxy powder coat finish)



Asepsis 21[®] Doctor's Cart

- Three automatic Kink Valve handpiece positions
- 9 5/8" x 13 5/8" stainless steel instrument tray
- Adjustable height "H" frame (epoxy powder coat finish)



Asepsis 21[®] Hygienist's Cart

- Two automatic Kink Valve handpiece positions
- Assistant's instrumentation: HVE, saliva ejector and syringe
- 9 5/8" x 13 5/8" stainless steel instrument tray
- 13" x 18" optional work surface (shown here)
- Adjustable height "H" frame (epoxy powder coat finish)