



INTEGRATED ACCESSORIES



Sopro Intraoral Camera

- Sopro Touch allows for easy image capture without blurring
- Light weight and easy to clean
- Analog and digital (USB 2.0) compatible

617 Camera (shown above)

- Wide 105 degree angle of view (easy to see distal areas)
- Curved design which follows the natural shape of the mouth
- High fidelity images
- Latest generation of LEDs utilizing condensers to focus light with minimal shadowing
- Fixed-focus and large depth of field (portrait down to a single tooth)

717 Camera

- SoproShade allows for easier shade matching and information sharing
- Highly sophisticated optical unit provides unsurpassed image quality
- Large depth of field (infinite down to 1mm macro)
- New generation LED technology
- Thinnest distal end on the market



SATELEC® Mini LED Curing Light

- Second generation LED curing light combining power, efficiency and speed
- Only 6 to 12 seconds to polymerize 3mm of any composite
- Wide light spectrum (420-480 nm)
- Three modes of cure (fast, pulse and ramp)
- Comes standard with light shield and autoclavable 7.5 mm light guide
- Silent operation (no fans)
- Stylish handpiece is compact, light weight and well balanced
- Shock resistance and easy to clean



SATELEC® SP Newtron™ Ultrasonic

- 5th generation Piezo ultrasonic instrument
- Compact, light weight and autoclavable
- Powerful, yet gentle for applications beyond scaling such as perio, endo and surgery
- Patented "Cruise Control" for automatic regulation of frequency and power
- Color coding system to ensure proper tip and power setting
- Most extensive range of tips in the market, over 80
- Introductory kit includes 3 tips (#1, 10Z, H3)
- LED version is also available (shown above)



Dentsply Cavitron® SPS G98A

- The industry's most popular ultrasonic scaler
- Magnetostrictive technology
- SPS technology automatically adjusts system power to maintain tip stroke, compensating for load and conditions
- Water flow adjustment conveniently located on the handpiece
- 30 KHz insert included



Bien-Air Optima MX INT

- The #1 manufacturer of electric micromotors in the world
- Brushless, autoclavable motor
- Perfectly silent and very gentle to operate
- Motor speeds of 100 to 40,000 RPM for a broad range of procedures
- Programmable high/low RPM and low RPM torque limits with auto reverse
- Quick-connect coupling rotates over 360 degrees
- Accepts E type attachments
- Available in one or two motor systems
- Three year warranty

Midmark is an ISO 13485 and ISO 9001 Certified Company.
For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495
Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214
or visit our website at midmark.com

© 2009 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0804-00 Rev.A1 (9/09)