

**EyeMag Pro - Medical Loupes**  
**Seeing is believing.**



**EyeMag® Pro F and EyeMag® Pro S Medical Loupes**  
**Sophisticated optics for higher visualization standards.**

Optimal treatment requires optimal visualization. EyeMag Pro loupes from Carl Zeiss generate a brilliant, stereoscopic image facilitating better differentiation of ultra-fine structures. The result: increased efficiency and elevated levels of care. Select your loupe from a wide range of working distances and magnification levels. Carl Zeiss EyeMag Pro loupes offer a choice in carrier system: a titanium frame (F) or a headband system (S). You'll be amazed by the unparalleled optical quality of the EyeMag Pro system.



*EyeMag Pro F*

*EyeMag Pro S*

Optical systems					
Working distance (mm)	500	500	450	450	400
Working distance (inch)	20	20	18	18	16
Magnification	3.2x	4x	3.3x	4x	3.5x
Field of view (mm)	115	93	100	81	86
Working distance (mm)	400	350	350	300	300
Working distance (inch)	16	14	14	12	12
Magnification	4.3x	3.6x	4.5x	4x	5x
Field of view (mm)	68	71	56	56	44

Compliance for EyeMag Pro F and EyeMag Pro S Medical Loupes and accessories:



EyeMag is a registered trademark of Carl Zeiss.

**Carl Zeiss Surgical GmbH**

A Carl Zeiss Meditec Company  
 Carl-Zeiss-Straße 22  
 73447 Oberkochen  
 Germany  
 Fax: +49 (0) 7364/20-4823  
 Email: [surgical@meditec.zeiss.com](mailto:surgical@meditec.zeiss.com)  
[www.meditec.zeiss.com/loupes](http://www.meditec.zeiss.com/loupes)  
[www.meditec.zeiss.com/contacts](http://www.meditec.zeiss.com/contacts)

**Quality optics and sophisticated design**

**Precise vision ...**

- High image quality and color fidelity – brilliant image that extends to the peripheral zones
- High depth of field yields excellent depth orientation
- Large fields of view allow for an optimum overview of the working field

**Work in comfort ...**

- The ergonomic design of the carrier system ensures an outstanding fit
- Customizable settings promote an upright posture regardless of the viewing angle
- The large selection of working distances permits an optimal treatment position

**Intuitive operation ...**

- Simple adjustments for fast, easy operation
- Flip-up function facilitates unobstructed vision and easy communication with the patient
- Sterilizable contact guards allow quick flip up/down operation even under sterile conditions

**EyeMag Pro (Kepler type) Medical Loupes include:**

- EyeMag Pro F telescopes and carrier system (titanium eyeglass frame available in three sizes)
- EyeMag Pro S telescopes and carrier system (adjustable headband)
- Protective lens caps
- Sterilizable contact guards
- Optics cleaning cloth

