

ELEVANCE  
CHAIR



E L E V A T I N G  
P E R F O R M A N C E

The Elevance™ Dental Chair – sophisticated styling, optimal patient access and exceptional patient comfort.

The result of years of research and development, this patient chair is a genuine departure from any chair on the market today.

The unique Cantilever Forward™ design accentuates the dramatic styling, allows for unsurpassed range of travel and offers greater flexibility to enhance the way you work, whether seated or standing. Though elegant, the Elevance™ chair is supported by a powerful hydraulic system and cast aluminum backrest to support patients weighing up to 450 lbs.

The thin, narrow backrest provides optimal access to the patient. The sculpted backrest and exclusive integrated armrests support the upper body while the patient enjoys a relaxing massage from the heat and massage system, reducing patient anxiety.

For those who demand the very best...elevate your expectations, choose Elevance™.



## Greater Access Without Compromise

Optimal access to the oral cavity is crucial for maintaining an ergonomically correct working posture and maximizing performance. The Elevance™ chair has a unique backrest that is thin and narrow and cradles patients with exclusive integrated armrests for exceptional patient comfort and access. Providing an unprecedented range of travel, the Cantilever Forward™ design brings the patient to you and is conveniently controlled with the industry's first wireless chair control.



Starting at a low seat height of 15" and extending to a seat height of 34", the Cantilever Forward™ design delivers the greatest range of travel of any chair available today. While seated or standing, the Elevance™ chair offers unsurpassed flexibility to enhance the way you work.



### Magnetic Headrest

The optional, double articulating, magnetic headrest comfortably supports and accurately positions the patient for full access.




### Wireless Chair Control

Offering you greater freedom, the industry's first wireless chair control includes four programmable presets and a composite safe mode for the dental light. (optional)



### Rotation Brake Release

Maximize your operatory space with convenient access to electronic chair rotation release pedals, intuitively located in the chair base.



Concave footrest nestles feet together, promoting relaxation

Contoured seat pocket for proper positioning and increased comfort

Advanced lower lumbar support for the lengthiest of procedures

Fully integrated heat and massage (optional)

High density foam for increased support and long life

Optional Ultraleather upholstery for luxurious softness

Integrated armrests cradle the upper body and allow for ease of entry/exit

Cast aluminum backrest for greater support and patient confidence

Sculpted headrest cradles the patient for greater comfort and stability

## Exceptional Comfort

Patient comfort is one of the most important factors in creating an ideal patient experience. Building from Midmark's reputation for unsurpassed patient comfort, we've refined the patented, integrated armrest design which uniquely cradles the upper body for exceptional comfort. Further enhancing their experience, patients can now enjoy a relaxing massage from the new heat and massage system.



Synchronizing the movement of the backrest to the human body improves comfort, eliminates the need for the patient to readjust and allows the dental team to begin the procedure more quickly.

S E R E N I T Y



The SerenEscape™ system allows you to enhance the patient experience. Patients are offered a serene escape with a relaxing massage, while lumbar heat warms muscles for further relaxation. Personal preferences are easily adjusted with the patient's remote control.



User-friendly remote allows patients to adjust settings to their preference. Settings include high and low heat, three different massage patterns, adjustable speed and intensity. The system also includes an auto shut-off mode.

Elevance™ Procenter LR Operator



Elevance™ Procenter Console Operator



### Elevance™ Procenter Operator

The Elevance™ Procenter Operator embodies the perfect balance of patient access and comfort paired with reliable instrument delivery and premium dental lighting. Sophisticated styling and uncompromising quality create superior value.



### Ergonomically Refined

Midmark operator stools deliver premium styling and exceptional comfort in an ergonomic seating system. Multiple seat sizes, seat heights and unsurpassed adjustability allow you to better personalize your seating solution. The seating system offers a contoured seat, forward tilt and waterfall front edges, promoting a more balanced posture and improved comfort.

We have tried our very best to accurately represent product colors in this brochure. However, due to the nature of printed media, the color of actual products may differ slightly from how they appear on your particular printed piece. In order to confirm your color selection, we strongly recommend that you contact Midmark Customer Service at 1-800-MIDMARK to request a sample before placing your order.

**Midmark is an ISO 13485 and ISO 9001 Certified Company.**  
For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495  
Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214  
or visit our website at [midmark.com](http://midmark.com)

© 2009 Midmark Corporation  
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286  
Products subject to improvement changes without notice  
Litho in U.S.A. 007-0814-00 Rev. A1 (9/09)



Because we care.