

ULTRACOMFORT®  
ULTRATRIM®



## C O M F O R T

The Midmark UltraComfort® Dental Chair offers unparalleled patient comfort. With its fully contoured back and specially designed lumbar support, the UltraComfort® Chair will keep patients supported and comfortable through even the lengthiest of procedures.

Available in your choice of seamless Silhouette, plush Ultraleather and a limited seams version of Ultraleather upholstery.





Midmark's renowned quality starts at the core with superior components. A cast iron lift structure provides outstanding strength and stability. The robust, cast aluminum backrest cradles the patient, delivering a comfortable, secure patient, and the hydraulic drive system provides smooth, quiet, precise operation.

## Features

The patented, integrated armrests not only provide essential support, but allow for unencumbered patient entry and exit. The double articulating headrest comfortably supports and accurately positions the patient for full access while the optional magnetic headrest provides additional neck support and permits side to side movement of the head. Multiple chair controls are standard, featuring four programmable presets that make placing the chair and patient in your preferred position a cinch.

## A C C E S S

The Midmark UltraTrim® Dental Chair offers exceptional clinician access. We carefully plotted essential support areas for a wide range of body types, which allowed us to remove materials in just the right places. The result is a chair back that is narrower and 40% thinner than our standard chair back.

Available in your choice of seamless Silhouette, plush Ultraleather and a limited seams version of Ultraleather upholstery.



Chair back design is key to proper clinician positioning and patient access. While the UltraComfort® chair back is thin and narrow, allowing you to work in an ergonomically correct position, the UltraTrim® chair back takes this to a new level with a back so thin and so narrow that even direct patient access is easily achieved.

We have tried our very best to accurately represent product colors in this brochure. However, due to the nature of printed media, the color of actual products may differ slightly from how they appear on your particular printed piece. In order to confirm your color selection, we strongly recommend that you contact Midmark Customer Service at 1-800-MIDMARK to request a sample before placing your order.

**Midmark is an ISO 13485 and ISO 9001 Certified Company.**  
For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495  
Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214  
or visit our website at [midmark.com](http://midmark.com)

© 2009 Midmark Corporation  
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286  
Products subject to improvement changes without notice  
Litho in U.S.A. 007-0816-00 Rev. B1 (11/10)



Because we care.