

Seating

Midmark Veterinary
Products and Services



What you need most in seating is comfort and support. Our seating products offer wide, amply padded seats for greater comfort while a strong base structure provides superb stability.



Let our seating products provide the comfort and support you deserve.

Midmark veterinary seating products were designed for unsurpassed durability, even under the most rigorous clinical conditions. With many choices, colors and accessories, we are certain to have a stool to meet your specific needs.



Midmark 277 Air Lift Stool



Midmark 279 Air Lift Surgeon Stool

Classic Series Stools feature a polished aluminum/black accented 5-leg, star-base with heavy duty hooded casters, up to 350 lb load capacity and a variety of backrest and lift options.



276

Midmark Dental Stools are not only engineered to deliver premium styling and exceptional comfort, they are also ergonomically refined to address musculoskeletal pain common when performing procedures. Promoting a more balanced posture, the seating system offers a contoured seat with raised front to improve weight distribution and overall comfort. Comfort is also enhanced with waterfall front edges that reduce pressure on the back of the thighs. Multiple seat sizes and seat heights, paired with unsurpassed adjustability, allow you to better personalize your seating solution.

Ergonomically Refined

- a. Contoured seat with raised front to improve weight distribution and overall comfort
- b. Waterfall front edges reduce pressure on the back of the thighs
- c. One lever adjustment for locking seat and backrest tilt
- d. Ratcheted backrest for intuitive height adjustment without knobs
- e. Personalize your stool by selecting from multiple stool sizes, cylinder heights 15.25" - 27.5" and upholstery levels and colors
- f. Integrated footrest that automatically adjusts with seat height



Value Series Stools feature a black composite 5-leg star-shaped base with heavy duty hooded casters, up to 300 lb. load capacity and are available with backrest and manual or automatic adjustment from 18" to 24.5" in height.






Optional armrests for the **Dental Doctor's stool** offer support for elbows in the neutral shoulder position to reduce fatigue of the neck and shoulder muscles.



Arm pad adjustments include height, swivel and slide for optimal positioning of the armrests to better accommodate your personal needs.



The **Dental Technician's stool** offers a ratchet release body support with an integrated footrest that automatically adjust with the seat height.



The 680 Side Chair comes with all the extras; contoured seat and back section for added support, sturdy, powder coated black steel tube frame with leveling feet for increased stability and a "wall-saver" kick back leg design that allows the chair to be close to the wall without damaging it. Available with or without high-impact plastic arms.



- Wall saver, kickback rear leg design
- Leveling feet
- 400 pound weight capacity



Classic Series Stool Specifications

Seat height:
274 - 278
18" - 24.5" (45.7 - 62.3 cm)

279
19" - 25.5" (48.3 - 64.8 cm)

Seat diameter:
16.25" (41.3 cm)

Seat cushion thickness:
3.25" (8.3 cm)

Back cushion dimensions:
275, 277
14" W x 7.5" H (35.6 x 19.1 cm)

Caster base diameter:
21.5" (54.6 cm)

Weight capacity:
274, 276
350 lbs. (158.8 kg)
275, 277 - 279
225 lbs. (115.7 kg)

Value Series Stool Specifications

Seat height:
270 - 273
18" - 24.5" (45.7 - 62.2 cm)

Seat diameter:
16.25" (41.3 cm)

Seat cushion thickness:
3.25" (8.3 cm)

Back cushion dimensions:
271 - 273
14" W x 7.5" H (35.6 x 19.1 cm)

Caster base diameter:
23" (58.4 cm)

Weight capacity:
270, 272
300 lbs. (136.1 kg)
271, 273
225 lbs. (102 kg)

Classic and Value Series Options

Locking casters
Soft rubber casters
Glides
6" Height extension
Foot ring

680 Side Chair Specifications

Seat height:
17" (43.2 cm)
Seat cushion dimensions:
18" W x 19.5" D
(45.7 x 49.5 cm)
Seat cushion thickness:
2" (5.1 cm)
Back cushion dimensions:
18" W x 13.5" H (45.7 x 34.3 cm)
Inside arm width: 17.5" (44.5 cm)
on model with arms
Outside arm width: 23.5"
(59.7 cm)
on model with arms
Weight capacity:
400 lbs. (181.4 kg)

Doctor and Technician Stool Specifications








Seat size:
17.5" x 16.5" (44.5 x 41.9 cm)
standard
16.25" x 16" (41.3 x 40.6 cm)
small
Cylinder height:
16.5" - 22" (41.9 x 55.9 cm)
standard
15.25" - 18.3" (38.7 x 46.5 cm)
short
Short technician's stool cylinder height is not available with a footrest.

Available Colors








We have tried our very best to accurately represent product colors in this brochure. However, due to the nature of printed media, the color of actual products may slightly differ from how they appear on your particular piece.

To request color swatches, please contact the Midmark Veterinary Customer Service Department at 1-800-MIDMARK. *There will be an additional charge for Midmark Ultra-Premium.*







Midmark Premium Solid Upholstery Selections

 Blueberry
 Cashmere Blue
 Silver Sage
 Terra Cotta
 Fossil Grey
 Tea Green
 Black









Midmark Ultra-Premium Upholstery Selections

 Ultra-Pacific Blue
 Ultra-Sky
 Ultra-Spruce
 Ultra-Canyon
 Ultra-Arctic
 Ultra-Olivine
 Ultra-Raven Wing

Midmark Premium Pattern Upholstery Selections

 Belagio
 Capri
 Siena
 Treviso
 Firenze
 Tuscany

Ritter Upholstery Selections

 Navy
 Dusty Blue
 Moss
 Clay
 Pebble Grey
 Shadow
 Perfect Plum
 Black

Midmark is an ISO 13485 and ISO 9001 Certified Company.
For more information or a demonstration, contact your Midmark dealer or Call: 1-800-MIDMARK Fax: 1-877-725-6495
Outside the U.S.A. Call: 1-937-526-3662 Fax: 1-937-526-8214
or visit our website at midmark.com

© 2009 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0778-00 Rev. A1 (8/09)



Because we care.